



Foetal Bovine Serum

Manufactured by the only Australian owned company supplying FBS in Australia



100% Australian Origin FBS



Completely traceable and independently quality tested



Manufactured in a fully validated State-of-the art facility



Certified by the European Directorate of Quality Medicines



Low Endotoxin and Haemoglobin Levels



ITEM CODE

SFBS

DESCRIPTION

Sterile, Filtered Foetal Bovine Serum

VOLUME

500ml

Now exclusively available from:

Interpath Services
Leading the way

ABN 79 007 196 581

Contact: **Customer Service**
Melbourne Office: **03 9457 6277**
Freecall: **1800 62 63 69**
E-mail: **sales@interpath.com.au**

TECHNICAL INFORMATION AND SPECIFICATIONS

QUALITY ASSURANCE:

FBS testing is completed by external, accredited laboratories to provide an independent guarantee in relation to the reporting of finished product test results.

Foetal Bovine Serum is tested for:

- pH
- Haemoglobin
- Total IgG
- Osmolality
- Total protein
- Biochemistry
- Endotoxin Levels
- Sterility
- Mycoplasma
- Virus isolation



Reports are available for customers to verify product quality and only those lots that meet strict requirements are released for sale.

Foetal bovine Serum is filtered through a series of sterile filters, starting at 5µm and finishing with three 0.1µm filters in a fully validated and dedicated state-of-the-art sterile filtration facility operated to cGMP standards.

ORIGIN and TRACEABILITY:

Australian derived bovine biologicals do not carry the risk of serious animal diseases such as Bovine Spongiform Encephalopathy (Mad Cow Disease), Foot and Mouth Disease and Bluetongue that can be found in North and South America, Europe and Asia.

Complete traceability of the origin of Bovogen Biological's FBS is available from farm to serum via the mandatory animal identification system used in Australia.



Bovogen Biologicals supply a 100% Australian made and owned product. As they are a member of the Australian Made campaign they are authorised to use this logo.

For a free sample of Foetal Bovine Serum from Interpath Services,
please complete the attached information and return to your local
representative, or fax to:

Melb: (03) 9458 4010

Freefax: 1800 658 210

Interpath Purchase Order No.: _____

Institute Name:.....

Department:.....

Contact Name:.....

Delivery Address:.....

.....

Telephone Number:.....

E-Mail Address:.....

Estimated total quantity of bottles required:.....

Number of samples required:.....

OFFICE USE ONLY:

Batch No:.....

Date Dispatched:.....

ITEM CODE

SFBS

DESCRIPTION

Sterile, Filtered Foetal Bovine Serum

VOLUME

500ml

Interpath Services
Leading the way

ABN 79 007 196 581

Contact: **Customer Service**

Melbourne Office: **03 9457 6277**

Freecall: **1800 62 63 69**

E-mail: **sales@interpath.com.au**