

# Freloads®



A simple, unique, cost effective way to purchase and renew your racked pipette tips.

Ask for Freloads® by name and say goodbye to costly, space wasting racks forever!

**Free** up your bench space

**Free** up having to store and dispose of empty racks

**Free** your funds for other purchases

**Free** yourself from dubious quality concerns







**Free** up landfill waste

Unit 1 46 Sheehan Road  
West Heidelberg VIC 3081  
Ph: 03 9457 6277  
Fax: 03 9458 4010

Freecall: 1800 62 63 69  
Freefax: 1800 65 82 10  
E-Mail: [sales@interpath.com.au](mailto:sales@interpath.com.au)  
[www.interpath.com.au](http://www.interpath.com.au)

**Interpath Services®**  
Leading the way

ABN 79 007 196 581

		COMPATIBLE WITH							
SIZE	PART NUMBER		Socorex	Biohit	Brand Transferpette	Eppendorf	Finnpipette	Gilson Pipetman	Nichiryo/Oxford
<b>10</b>	<b>33400</b>	Low Retention Freload® Tips	Acura 8/12 Channel 0.5 - 10µl Acura 0.5 - 10µl	Electronic 0.2 - 10µl Proline 0.5 - 10µl	0.1 - 10µl	Electronic 0.5 - 10µl Reference 0.1 - 10µl 0.5 - 10µl Research 0.5 - 10µl	0.2 - 10µl	P2 P10	Benchmate 0.5 - 10µl
	<b>39700</b>	Low Retention Bulk Tips							
 <p>Tip length: 31.25mm</p>									
<b>10 XL</b>	<b>33410</b>	Low Retention Freload® Tips	Acura 0.1 - 2µl 0.5 - 10µl Acura 8/12 Channel 0.5 - 10µl	Electronic 0.2 - 10µl Proline 0.5 - 10µl	0.5 - 10µl	Reference 0.1 - 2.5µl 0.5 - 10µl Research 0.1 - 2.5µl 0.5 - 10µl	0.5 - 10µl	P2 P10	Benchmate 0.1 - 2µl 0.5 - 10µl
	<b>39710</b>	Low Retention Bulk Tips							
 <p>Tip length: 45.84mm</p>									
<b>200</b>	<b>33430</b>	Low Retention Freload® Tips	Acura 20 - 200µl Acura 8/12 Channel 20 - 200µl 40 - 350µl	Electronic 5 - 100µl Proline 5 - 50µl 50 - 200µl	10 - 100µl 20 - 200µl	Electronic 5 - 100µl 2 - 200µl Reference 10 - 100µl 20 - 200µl Research 10 - 100µl 20 - 200µl	Single 40 - 200µl Multi-channel 50 - 300µl	P200	Benchmate 10 - 100µl 20 - 200µl
	<b>39730</b>	Low Retention Bulk Tips							
 <p>Tip length: 59.54mm</p>									
<b>300</b>	<b>33450</b>	Low Retention Freload® Tips	Acura 20 - 200µl	Electronic 5 - 100µl 20 - 200µl Proline 10 - 100µl 20 - 200µl	20 - 200µl 15 - 300µl	Reference 50 - 200µl Research 20 - 200µl 30 - 300µl	Single 5 - 50µl 20 - 200µl 30 - 300µl Multi-channel 5 - 50µl 30 - 300µl	P200	Benchmate 20 - 200µl 50 - 300µl
	<b>39750</b>	Low Retention Bulk Tips							
 <p>Tip length: 59.44mm</p>									
<b>1000</b>	<b>33470</b>	Low Retention Freload® Tips	Acura 100 - 1000µl	Proline 200 - 1000µl	100 - 1000µl	Electronic 50 - 1000µl Reference 100 - 1000µl Research 100 - 1000µl	Single 100 - 1000µl 200 - 1000µl	P1000	Benchmate 100 - 1000µl
	<b>39770</b>	Low Retention Bulk Tips							
 <p>Tip length: 83.94mm</p>									
<b>1250</b>	<b>33490</b>	Low Retention Freload® Tips	Acura 100 - 1000µl	Proline 200 - 1000µl	100 - 1000µl	Reference 100 - 1000µl Research 100 - 1000µl	Single 100 - 1000µl 200 - 1000µl Multi-channel 100 - 1200µl	P1000	Benchmate 100 - 1000µl
	<b>39790</b>	Low Retention Bulk Tips							
 <p>Tip length: 98.30mm</p>									
Tip photographs not actual size.									

All Freload® tips up to 1000µl are packed in 10 trays of 96 tips. The 1000µl and 1250µl Freload® tips are packed in 8 trays of 96 tips. All bulk tips are packed in bags of 1,000.

#### RACKS

- 33497** Small tip rack to fit 10µl and 10XL Freload® tips
- 33498** Medium tip rack to fit 200µl and 300µl Freload® tips
- 33499** Large tip rack to fit 1000µl and 1250µl Freload® tips



— Certified —  
RNase/DNase  
— FREE —