

JENCONS

# Powerpette® Pro Pipette Controller

NEW

- › Fills a 25ml pipette in under 3 seconds on its fastest setting
- › Rechargeable during use
- › Powerful pump
- › 'Low battery' light
- › Environmentally friendly nickel Metal hydride battery
- › Suitable for all types of pipettes 1 - 100ml
- › Unique stand suitable for bench or wall mounting
- › CE marked
- › 1 Year Warranty



**Interpath Services**

Leading the way



Please Contact:  
Customer Service

Phone: 03 9457 6277

Freecall: 1800 62 63 69

E-mail: [sales@interpath.com.au](mailto:sales@interpath.com.au)

ABN: 79 007 196 581

## Powerpette® Pro Pipette Controllers

The latest addition to the Powerpette® family of pipette controllers, the Powerpette® Pro fits the hand comfortably with the pipette holder at an angle ideal for use within a laminar flow or microbiological safety cabinet.

The Powerpette® Pro Pipette controller delivers efficient performance with a powerful yet quiet motor to speed up large volume pipetting.

The Powerpette's mode selection switch enables selection of high or low modes both allowing variable aspirate and dispense speeds (with blow out). The gravity dispense mode is designed for use with 'To Deliver' (TD) pipettes. The speed of suction and dispensing is controlled through the concave finger triggers, designed to provide a comfortable and positive grip and to require the minimum of effort.

The Powerpette® Pro is supplied with 0.45µm and 0.2µm hydrophobic filters, as standard, for effective protection against airborne contamination of the unit during aspiration, and of your samples during dispensing.

### Specifications

Batteries	2 x 1.2V nickel metal hydride
Battery service life	4 hours continuous life
Charger	75mA/h for 16 hours
Overall weight	180g without stand/wall bracket
Filter, µm	Hydrophobic 0.45 or 0.2
Case material	ASA
Nosecone material	ASA
Pipette holder	Silicone
Pipette types, ml	Plastic or Glass 1 - 100

Cat. No.	Description
266-185*	Powerpette® Pro Pipette Controller

The Powerpette® Pro is fitted with a 0.45µm filter and supplied complete with charger, two spare hydrophobic filters (1x0.45µm and 1x0.2µm) bench stand/wall bracket and instruction manual.

\* This catalogue number refers to units supplied with a Australian compliant charger.



### Accessories

Cat. No.	Description	Pack Size
266-003	Replacement filter Set 0.45µm	5
266-089	Replacement Filter Set 0.2µm	5
266-004	Replacement Silicone Pipette Holder	1
266-090	Sterile filter 0.2µm	1
266-186	Spare bench stand for Powerpette® Pro	1



The bench stands for Powerpette® Pro can be used to house the controller with or without a pipette in place. The controller can be charged in either position.