



## TECHNICAL DATA SHEET

### **PRODUCT DESCRIPTION**

<b>PRODUCT CODE</b>	155C
<b>GENERAL DESCRIPTION</b>	TAMPONE SECCO A/PL 1000 PKG
<b>DETAILED DESCRIPTION</b>	10 boxes each containing 100 tubes with white cap with inside a plastic applicator with rayon tip in bulk

### **PRODUCT COMPONENTS**

<b>APPLICATOR</b>	PLASTIC APPLICATOR	<b>DIMENSIONAL SPECIFICATIONS</b>
		STICK LENGTH: 144.0 mm $\pm$ 1.0 mm PLASTIC TIP DIAMETER: 2.50 mm $\pm$ 0.20 mm FIBRE TIP LENGTH: 15.0 mm $\pm$ 2.0 mm FIBRE TIP DIAMETER: 4.50 mm $\pm$ 0.50 mm TOTAL LENGTH: 147.0 mm $\pm$ 1.0 mm
<b>CAP</b>	WHITE CAP	<b>DIMENSIONAL SPECIFICATIONS</b>
		EXTERNAL CAP DIAMETER: 12.00 mm $\pm$ 0.10 mm CAP LENGTH: 32.80 mm $\pm$ 0.20 mm INSERTION DIAMETER: 9.30 mm $\pm$ 0.05 mm
<b>TRANSPORT TUBE</b>	TUBE	<b>DIMENSIONAL SPECIFICATIONS</b>
		TUBE LENGTH: 154.50 mm $\pm$ 1.00 mm INSERTION DIAMETER: 9.0 mm $\pm$ 0.2 mm EXTERNAL DIAMETER: 12.0 mm $\pm$ 0.5 mm

This bulletin (Rev.02-27/06/2014) is for your guidance and is based on information believed to be reliable and should be periodically reviewed.



## TECHNICAL DATA SHEET

### **PACKAGING SPECIFICATION**

<b>INNER BOX QTY</b>	100
<b>OUTER BOX QTY</b>	1000
<b>PACKAGING</b>	BULK

### **IMPORTANT INFO**

<b>STERILIZATION METHOD</b>	GAMMA IRRADIATION
<b>STORAGE TEMPERATURE</b>	5/25°C-DO NOT FREEZE
<b>SHELF LIFE</b>	36 MONTHS