



PRODUCT SPECIFICATION SHEET

Product Code: BZFMELB/2A
Product Description: Face Mask, Procedural, 3-Ply, Blue, Level 2
Brand: Premier

Product Information

Type: Procedural Face Mask, Level 2
Material: Nonwoven PP, Spandex
Colour: Blue
Design & Features: Blue, 3-Ply, Pleated, Ear-Loop, Latex-Free
Barrier Protection Level: (AS) 4381: 2015 - Level 2
Bacterial Filtration Efficiency (BFE) %: $\geq 98\%$
Differential Pressure (Delta P) mm H₂O/cm²: < 5.0
Synthetic blood penetration resistance: 120mm Hg
(min pressure in mm Hg for pass result)
Intended Use: Healthcare industry, food industry, general cleaning, single use only
Quality Assurance: ARTG Number - 333127
ISO 13485:2012, CE
Shelf Life: 3 years from date of manufacture
Storage: Store in a cool, dry area, away from direct sunlight and humidity
Country of Manufacture: China

Packing Information

UOM: 50 pieces/box, 40 boxes/carton
Inner Box Dimensions (mm): 180 x 100 x 80
Outer Carton Dimensions (mm): 520 x 380 x 340
Carton Weight (kg): 6.5
Inner Barcode: 9325191032724
Outer Barcode: 19325191032721

Image:

