

WASTE BAGS



BIOHAZARD

FOR AUTOCLAVING BIOHAZARD WASTE. AVAILABLE PRINTED WITH BIOHAZARD LABEL OR UNPRINTED



Item No.	Size	Thickness	Heat	Colour	Ctn Qty
ACB001	27 x 63cm	50µm	121°C	Clear with Print	500
ACB002	55 x 70cm	50µm	121°C	Clear with Print	250
ACB003	72 x 91cm	50µm	121°C	Clear with Print	250
ACB051	27 x 63cm	50µm	121°C	NO Print / Clear	500
ACB052	55 x 70cm	50µm	121°C	NO Print / Clear	250
ACB082	55 x 70cm	80µm	131°C	Clear with Print	250
ACB083	72 x 91cm	80µm	131°C	Clear with Print	100
ACB091	27 x 63cm	50µm	121°C	Yellow with Print	500
ACB092	55 x 70cm	50µm	121°C	Yellow with Print	250
ACB093	72 x 91cm	50µm	121°C	Yellow with Print	250
ACB098	120 x 85cm	50µm	121°C	NO Print / Yellow	100

Yellow Clinical Waste Bags

Item No.	Description	Ctn Qty
YIW5060	500mm x 600mm	500
YIW6090	600mm x 900mm	250
YIW63585	635mm x 850mm 50µm	250
YIW6691	660mm x 910mm 75µm	200
YIW7010	700mm x 1000mm	250
YIW130LT	860mm x 1220mm	200
YIW240LT	1130mm x 1450mm	100




Purple Cytotoxic Waste Bags

Item No.	Description	Ctn Qty
CYT5060	500mm x 600mm	100
CYT6380	630mm x 800mm	200
CYT7190	710mm x 900mm	100
CYT140LT	600 + 400mm x 1450mm (Gusset)	100
CYT240LT	590 + 590mm x 1500mm (Gusset)	100/Roll



SPECIMEN BAGS

Medication Transfer Bag



Item No.	Description	Ctn Qty
MS050	Green Ziplock Printed Patient Own Medication Transfer Bag	500
MS070	Red Ziplock Printed Ward Transfer Inpatient Medication Bag	500

Specimen Carrier Bags

Item No.	Description	Ctn Qty
SCPB165	Clear Specimen Carrier Bag	1000
SCPB166	Red Specimen Carrier Bag	1000
BSCP168	Blue Specimen Carrier Bag	1000



Ziplock Press Seal Bags Also Available

24 Sizes ranging from 50mm x 75mm up to 530mm x 550mm

Please contact your Sales Representative for further information and pricing.