


PDS No. 633102/ 633161/ 63318x	<b>PRODUCT DATA SHEET</b>				Page 1 of 1
Revision 10	<b>Petri Dish with Vents</b>				 <b>greiner bio-one</b>
	Greiner Item-No. 633102 / 633161 / 63318x				
Valid for Item-No.:	<b>633102</b>	<b>633161</b> (sterile)	<b>633180</b>	<b>633185</b>	

1.	Description / Specification	
1.1	Description	Petri dish with vents 633102 (94 x 16 mm): heavy design 633161 (94 x 16 mm): heavy design, sterile 633180 (94 x 16 mm): standard 633185 (94 x 15 mm): standard
1.2	Dimensions	See Customer Drawing 633102, -161: total weight: 17,0 – 18,5 g 633180: total weight: 13,0 – 14,7 g 633185: total weight: 13,0 – 14,5 g
1.3	Volume	Max. volume: 80 ml
1.4	Material / Resin	Dish: PS (Polystyrene), free of heavy metal Lid: PS (Polystyrene), free of heavy metal
1.5	Colour	Dish: clear Lid: clear
1.6	Sterilisation	633102, -180, -185: no 633161: SAL 10 <sup>-3</sup>
1.7	Quality Control	- <u>Raw Material-Control</u> : physical testing - <u>Product-Control</u> : testing of attributive and variable characteristics in accordance with the valid specification
1.8	Other Information	- For single use only - Elevations for ventilation of culture

2.	Features	
2.1	Basic features	Hydrophobic
2.2	Temperature range	-20°C to +60°C
2.3	Autoclavability	No
2.4	Centrifugation, max. RCF	N/A
2.5	Chemical Resistance	See homepage: <a href="https://www.gbo.com/en_INT/know-how-services/download-center.html">https://www.gbo.com/en_INT/know-how-services/download-center.html</a>
2.6	Shelf life	633102, -180, -185: n/a 633161: 5 years after month of production
2.7	Other Information	-

3.	Packaging	
3.1	Pieces / Bag	20
3.2	Pieces / Box	480
3.3	Lot-No.	F YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)
3.4	Other Information	-

4.	Other Information
	-

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	<b>CONFIDENTIAL:</b> Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 09	Date 2 March 2015	Date 3 March 2015	Date 3 March 2015	
Date 10.07.2013	Name S. Kaelberer	Name Dr. T. Schreiber	Name A. Schulz	

DISCLAIMER: The description of a certain product can only be considered as a guidance, because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.